# To Study the Effectiveness Of Different Common Oils (Castor Oil, Cotton Seed Oil, Coconut Oil, Kerosene Oil, Mustard Oil) In Forming Emulsions

# **Theory**

Depending upon the physical properties and the chemical composition of the oil, some oils form emulsions readily with water whereas many others form with great difficulty.

# **Apparatus**

Five stoppered bottles or boiling tubes, measuring cylinder, stop-watch or simple watch, 5 ml pipettes.

### **Materials Required**

Castor oil, cotton seed oil, coconut oil, kerosene oil, mustard oil and 1% soap solution or sodium oleate solution. (1 gm dissolved in 100 ml of distilled water).

#### **Procedure**

- 1. Take five stoppered bottles and wash them with water. Label them as A, B, C, D and E.
- 2. Take 5 ml of each of the oils separately in tubes A, B, C, D and E.
- 3. Add 50 ml of distilled water to each tube.
- 4. Take bottle A, stopper it and shake vigorously for one minute. Then allow it to stand.
- 5. Note the time taken for the two layers to separate out.
- 6. Similarly, take bottles B, C, D and E and note the time taken for the separation of two layers in each case.
- 7. Now add 10 drops of 1% soap solution or 1% sodium oleate solution to each of the five bottles and find out the time taken for the two layers to separate. Record the observations.

## **Observations**

Bottle or Tube	Name of the oil	Time taken for the separation of layers	
electron bank to c	u bojuga a sassa A s	Without soap/detergent	With soap/detergent
A	Castor oil	— s	— s
В	Coconut oil	— s	— s
c	Mustard oil	— s	— s
D	Cotton seed oil	— s	— s
Е	Kerosene oil	— s	<u>-s</u>

## Result

It is clear from the above observations that...... oil takes longest time to get separated from its emulsion and is rated 1 and...... oil takes the minimum time and is rated 5. The decreasing order of stability or effectiveness is.......

## **Precautions**

- 1. Add equal number of drops of soap solution to all the bottles.
- 2. Each bottle should be shaken vigorously and for same time. .
- 3. The time should be recorded carefully. Start the stopwatch immediately after shaking is stopped and stop it immediately when the two layers just separate.