

Exercise 10

Aim: To verify the Mendel's Law of Independent Assortment

Principle: In a dihybrid cross, the segregation of one gene pair is independent of the segregation of the other pair. It means that genes of two different traits assort independently to give a probability ratio equal to segregation probability ratio of one allele pair X segregation probability ratio of other allele pair, which comes to, $(3:1) \times (3:1) = 9:3:3:1$

Requirement: Plastic beakers; 64 plastic beads each of yellow, green, red and white to represent, yellow and green colour of seed coat and red and white flowers respectively and napkin/hand towel

Procedure

Students are to work in pair.

The following steps are to be followed sequentially:

- (i) Place 64 beads of each colour in four separate beakers.
- (ii) Put the beakers containing the yellow and red beads on your left side, and those containing the green and white beads on your right side. The beakers on your left side represent plants bearing yellow seed and red flower (dominant character YY, RR). Beakers on the right side represent plants bearing green seeds and white flowers (recessive character yy, rr). These are the two parental types having contrasting forms of two different characters.
- (iii) Stir the beads in each beaker with a pencil/pen. Each bead now represents alleles in the male and female gametes.
- (iv) Pick up one yellow, one green, one red and one white bead, and put them together on the napkin spread on the table.
- (v) Continue picking up and putting together of the beads of all colours as mentioned in the previous step, till all the beads are utilised.
- (vi) Note that in all, 64 such 4-bead clusters are obtained representing the F_1 individuals. Ascertain their genotype and phenotype.
- (vii) Next step is to cross these F_1 individuals to raise the F_2 generation. Let us suppose half of the 4-bead clusters (32 clusters) represent the male parents and the remaining half (32 clusters) the female parents. Now put the 32 red and 32 white beads together in one beaker (numbered-I), and similarly put 32 yellow and 32 green beads together in other

beaker (numbered-II). These two beakers represent F_1 female. Similarly put remaining 32 red + 32 white beads in beaker numbered-III, and 32 yellow and 32 green one in beaker numbered-IV to represent the F_1 male. The arrangement can be presented as below

Female F_1	Male F_1
32 red + 32 white (Beaker I)	32 red+32 white (Beaker III)
32 yellow + 32 green (Beaker II)	32 yellow + 32 green (Beaker IV)

- (viii) Stir the beads in each beaker with a pencil. In order to raise the F_2 generation, pick up (with eyes closed) one bead from the beaker-I of female and one bead from the beaker-III of the male, and put into the palm of the partner student. Similarly, pick up one bead each from the beaker-II of female and beaker IV of male to put in the palm of the partner. This partner would now keep all the four beads together (to represent the F_2 individual). Continue this process till all beads are utilised. At the end, 64 F_2 individuals (each represented by a 4-bead cluster) are obtained.
- (ix) Determine the genotype and phenotype of each of the 64 F_2 individuals and write down the number of individuals of different genotypes and phenotypes in the tabular form (given below), remembering that Y (yellow seed colour) is dominant over y (green seed) and R (red flower) is dominant over r (white flower).
- (x) Repeat the whole procedure (steps i to ix) six times, and tabulate your results.

Observation

Tabulate the results as follow:

Symbol (-) indicates the presence of corresponding dominant or recessive allele e.g. Y or y and R or r.

Summarise your results (adding together the data of all the six repeats)

F1 Generation

- (a) Total number of individuals: _____
- (b) Phenotype (s) _____
- (c) Genotype (s) _____

EXERCISE 10

Generation & repeat No.	Total No. of offsprings	Genotype				Phenotype			
		Y-R-	Y-rr	yyR-	yyrr	Yellow Red	Yellow White	Green Red	Green white
F 1									
1.									
2.									
3.									
4.									
5.									
6.									
Total									
F 2									
1.									
2.									
3.									
4.									
5.									
6.									
Total									

F2 Generation

- (a) Total number of individuals _____
- (b) Phenotypes _____
- (c) Number of individuals in each phenotypic class:
- | <i>Number</i> | <i>Phenotype</i> |
|---------------|------------------|
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
- (d) Phenotypic ratio _____
- (e) Genotypic ratio _____

- (f) Number of individuals of each genotypic class:

Number**Genotype**

- (g) Genotypic Ratio _____

Discussion

The four phenotypic classes in the F₂ generation are in ratio of 9:3:3:1 as expected from the **Law of Independent Assortment**. The genotypic ratio would be (1:2:4:2): (2:1):(2:1):1.

Note

1. In case six repeats of the experimental procedure are not feasible due to time limitations, either the number of repeats be slashed down to three or the data from single repeat of six different pair of students may be pooled together to make the final calculations.
2. This Law of Independent Assortment was later found to be true only for traits present on two different homologous pair of chromosomes, that is, the two are not linked together. The linked traits do not assort independently, rather they are inherited together (linked) except when crossingover separates them.
3. It is quiet likely that you may not find your data exactly in the expected ratio, instead almost approximate to it. The statistical significance of this deviation from the exact expected ratio due to probality can be checked using chi-square (χ^2) test, about which you will study in higher classes.

Questions

1. Linked traits fail to assort independently. Explain.
2. How is independent assortment of alleles important from the point of view of variation?